

Table 1. Summary of reports on the number of MscL channels per *E. coli* cell.

| Number per Cell | Method | Reference |
|-----------------|-------------------|-----------|
| 50 | radiolabeling | [24] |
| 10–100 | electrophysiology | [11] |
| 4–5 | electrophysiology | [16] |
| 10–15 | electrophysiology | [25] |

doi:10.1371/journal.pone.0033077.t001

11. Blount P, Sukharev SI, Moe PC, Martinac B, Kung C (1999) Mechanosensitive channels of bacteria. *Methods Enzymol* 294: 458–482.

16. Stokes NR, Murray HD, Subramaniam C, Gourse RL, Louis P, et al. (2003) A role for mechanosensitive channels in survival of stationary phase: regulation of channel expression by RpoS. *Proc Natl Acad Sci U S A* 100: 15959–15964.

24. Hase CC, Minchin RF, Kloda A, Martinac B (1997) Cross-linking studies and membrane localization and assembly of radiolabelled large mechanosensitive ion channel (MscL) of *Escherichia coli*. *Biochem Biophys Res Commun* 232: 777–782.

25. Booth IR, Edwards MD, Murray E, Miller S (2005) The role of bacterial ion channels in cell physiology. In: Kubalski A, Martinac B, eds. *Bacterial Ion Channels and Their Eukaryotic Homologs*. Washington, DC: American Society for Microbiology. pp 291–312.